

## Endomethasone Powder, 14g



Behind Every Dentist

### Indications

- Permanent root canal filling and sealing in dental procedures, with systematic insertion of points

### Features & Benefits:

- Neither resorbable once hardened nor retractable over the years that means it does not break down or move
- Antiseptic and anti-inflammatory activity lasts for several hours after being placed in the canal
- Radiopaque and easy to insert in to the canal
- Endomethasone has all the various above given qualities , owing to its exipients and antiseptic active substances polyoxymethaylene and thymol iodide



**Composition :** (1) Hydrocortisone Acetate 1g (2) Polyoxymethylene 2.2 g (3) Exipients : thymol iodide , E 110 , Barium Sulphate , Zinc Oxide , Magnesium Stearate q.sad 100 g

### Directions for Use

- The treated tooth should be perfectly isolated (using a cotton roll or a rubber dam).
- After carefully shaping the canal, cleaning it with an antiseptic solution, rinsing it and drying it with paper points, prepare the canal paste, progressively pouring the powder into eugenol until a texture similar to petrolatum is obtained.
- Recommended Mixing Ratio is : 3-6 drops of Endomethasone Liquid in 2-Spoonfuls of Endomethasone powder. (Spoon provided in packing).
- Then fill the root canal using a paste-carrier, and position into the paste one or several points, which have previously been adapted to the canal shape