

Sodium Hypochlorite 3% solution 250ml

Solution for root-canal irrigation



SWISS QUALITY

PROPERTIES AND INDICATIONS

Sodium Hypochlorite 3% solution has a solvent action on pulp tissue and organic debris, and is used for irrigation of root canals.



CONTRA-INDICATIONS

Sodium Hypochlorite solution is caustic and not suitable for application to wounds or areas of infection in soft tissues.

COMPOSITION

Sodium Hypochlorite (13-14% active Chlorine)~ 25%, Stabilizer 0,08%, demineralized water ad 100%

UNIVERSAL  DENTAL (PVT) LTD.